

Kwik-Stop Syptic Powder With Benzocaine

Kwik-Stop is an aid to stop bleeding caused by clipping nails, docking tails, and trimming beaks and minor cuts. Pressure bandaging to be used in conjunction with product following tail docking.

Directions: Apply with moistened cotton tipped applicator to the cut nail or other superficial bleeding area using moderate constant pressure for five to ten seconds. Do not use in deep wounds or body cavities, or on burns.

Contains: Ferric subsulfate, aluminum chloride, diatomite, bentonite, copper sulfate, ammonium chloride, benzocaine

For external veterinary use only to control bleeding for dogs, cats & birds.

Keep out of reach of children.

Packaging: 14 grams

Manufacturer: ARC Laboratories